THE PATHCARE NEWS

16s rRNA PCR A new molecular test improves PathCare laboratories' capacity to detect and identify bacteria involved in 'difficult' infections

A 16s rRNA molecular test is now offered by the molecular laboratory in our reference laboratory at N1 City, Goodwood. The test is useful in the diagnosis of bacterial infection in cases where culture does not yield positive results due to fastidious growth requirements of the infecting agent, very low organism counts or when organisms are no longer viable due to exposure to antimicrobial agents. This technique can also be utilized to identify organisms cultured when conventional methods fail.

Q: What is 16s rRNA?

A: It is ribosomal RNA which is conserved in all bacteria, but has highly specific nucleic acid sequences unique to every species.

Q: How does the test work?

A: A PCR is done on the specimen in order to 'fish' for any bacterial RNA which may be present. These segments are then amplified and sequenced in order to obtain a genetic 'fingerprint', which can then be identified with the use of a web-based databank.

Q: On what specimens can this be performed?

A: Any specimen from a normally sterile site is suitable (for example heart valves, CSF, and aspirates from vitreous fluid and joints). Care must be taken to avoid contamination, as any bacterial RNA which may be present in the sample will be detected.

Q: Which organisms can be detected and identified using the 16s rRNA technique?

A: All bacteria, including tuberculosis and organisms with incomplete cell walls e.g. rickettsia and mycoplasmas. Viruses and fungi, however, cannot be detected using this technique.

Q: What about cases where the patient has already received antimicrobials?

A: Because this test does not rely on the viability of organisms, it can still be positive providing that the patient's immune system has not yet mopped up bacterial residues.

Q: What may influence the sensitivity of the test?

A: A control is built into the test to indicate the presence of inhibitors of the PCR reaction which may be present in the specimen. A very low number of organisms or a diluted sample (for instance peripheral blood) may also be below the limit of detection of this test.

Q: What may influence the specificity of the test?

A: Contamination of the specimen by skinflora or environmental organisms may lead to failure of the test, or problems with its interpretation, as organisms not relevant to the infective process may be detected. The detection of a pathogen associated with the specific infection increases the positive predictive value of the result.

Q: Would this test be useful in all kinds of infection?

A: This test performs best when infection is caused by a single type of bacteria. In cases of, for example, endocarditis and incompletely treated meningitis it has been shown to improve patient outcome and also to limit the unnecessary use of especially broad spectrum antimicrobial drugs. Mixed infections would not yield good results. Soft tissue and intra-abdominal infections are mostly polymicrobial in nature and this test will not be suitable in such cases.

Q: Can the 16s rRNA test also indicate antimicrobial sensitivity?

A: No. However, when dealing with a known bacterial infection better choices can be made in line with published data regarding antimicrobial susceptibilities of the identified organism.

Q: What is the turn-around time of this test?

A: 16s rRNA batches will be run twice a week and results should be available within 2 to 5 days.

Q: How do I request this test?

A: Because this is a highly specialized and low-volume test, it is not included on the menu of the standardised request form. It can be requested in the "other test" text box, and we also invite you to discuss cases where you consider this test to be of potential value to the management of your patient with the clinical microbiologist on call.

Compiled By: Dr Elizabeth Wasserman, Clinical Microbiologist, PathCare

Tel: 021 596 3400

For further reading: Clinical Microbiological Reviews 2004; 17(4): 840 - 862; Journal of Microbiological Methods 2009; 76: 217 - 225

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